



Region 1

Region 2

Region 3

Region 4

Region 5

Region 6

cellTRAY CT-2600 Slide Properties	
Size	1 " x 3" x 4mm
Material	Polycarbonate
Cell Lifetimes	> 5 days ⁱ
Well Diameter	1500 µm
Well Depth	3000 µm
# Regions	6
# Wells / Region	1
Total Wells	6
Single Well Volume	16.9 µL
Regional Channel Volume	7.1 µL
Total Regional Volume	24 µL
Total cellTRAY Volume	144 µL

ⁱ Using cellTRAY Fluidics System CT-2000F.